

Dynamics of EGFR Mutational Load in Urine and Plasma Correlates with


Abstract ID # 5529

Introduction

- Non-small cell lung cancer (NSCLC) is a heterogeneous and dynamic disease where testing for key mutations is essential.
- With the emergence of clonal resistance, obtaining serial biopsies to assist in the real-time treatment decision-making has proven challenging.
- Noninvasive urinary circulating tumor (ct) DNA-based liquid biopsy approach can be used to detect and track cancer driver mutations for rapid diagnosis and disease monitoring.
- Here, using a highly sensitive ctDNA mutation detection platform, we demonstrate detection of *EGFR* mutations in urine and plasma samples obtained from patients with NSCLC.

Clinical Study Design

- This is a prospective observational study of patients with non-squamous, tissueconfirmed EGFR, KRAS or ALK-mutant NSCLC preparing to receive a systemic treatment regimen.
- The primary endpoints are correlation between ctDNA and tumor-based molecular results, and measurable change in ctDNA with response by RECIST v1.1. Urine and blood specimens are collected from patients at baseline, on treatment,
- and at progression for ctDNA analyses.
- This study has enrolled 34 patients, and presented here is an interim analysis of 20 NSCLC patients with *EGFR*-positive tumors



*One patient with co-occurring EGFR and KRAS mutation.

Treatment Response in Advanced NSCLC

Jhanelle E. Gray¹, Ben Creelan¹, Tawee E. Tanvetyanon¹, Scott Antonia¹, Karen Johnson¹, Christine Sigua¹, Michelle M. Motschman¹, Jongphil Kim¹, Richie Reich¹, Braydon Schaible¹, Cecile-Rose T. Vibat², David Gustafson², Sandeep C. Pingle², Vlada Melnikova², Mark Erlander², and Eric B. Haura¹ ¹Moffitt Cancer Center, Tampa, FL, USA; ²Trovagene, Inc., San Diego, CA

N=34
66 (45-78)
15 (44.2%) 19 (56.1%)
6 (17.7%) 28 (82.5%)
26 4 4 0
4 30
4 15 15
25* 8* 2

EGFR Mutation Enrichment NGS Assay


- Mutant allele-enrichment PCR followed by NGS, utilizes a short footprint amplicon (42-46 bp) and selectively amplifies mutant EGFR DNA fragments.
- Quantitation of mutant DNA is achieved by including a set of standard samples in addition to clinical samples and controls to each single multiplexed NGS run.
- Mutation enrichment results in approximately 3000-fold increase in ratio of mutant over wild-type signal for low copy number inputs.
- Lower limit of detection: 1 copy in 18,181 WT GEq (0.006%) for EGFR Ex19del and L858R assays; and 2 copies in 18,181 GEq (0.01%) for EGFR T790M assay.

Correlation of ctDNA *EGFR* with Therapeutic Response


- 17 of 20 (85%) NSCLC patients with EGFR-positive tumors had detectable <u>concordant EGFR mutation</u> in pre-treatment (baseline) urine and/or plasma.
- 9 of 11 patients with matched serial ctDNA samples demonstrated detectable *EGFR* mutation signal in pre-treatment urine and/or plasma samples.
- Of 9 patients, 3 received single EGFR TKI, 3 received combination TKIs, 1 received chemotherapy, 1 received immune checkpoint inhibitors alone, and 1 received immune checkpoint inhibitors in combination with TKI.
- In all 9 patients, changes in ctDNA levels from baseline to every cycle on therapy were examined, and correlated with best response by RECIST v1.1.


Longitudinal Monitoring of ctDNA EGFR in Non-Responders

Representative case showing dynamics of urine and plasma ctDNA EGFR in **NSCLC** patient with progressive disease



Results





- and correlate with patient outcomes.

- For more information, please contact: Jhanelle Gray, MD 12902 Magnolia Drive, FOB1 Tampa, FL 33612 Jhanelle.Gray@moffitt.org



*Green arrow depicts change in ctDNA EGFR activating mutations from baseline (B) to cycle 2 day 1 (C2).

Summary

In both responders and non-responders, change in ctDNA EGFR activating mutations from baseline (B) to cycle 2 day 1 (C2, depicted by green arrow) correlated with best response by RECIST v1.1.

A 100% decrease in urine and plasma ctDNA EGFR activating mutations was observed in patients with partial response (PR, n=3) or stable disease (SD, n=3).

A <90% decrease or an increase in urine and plasma ctDNA EGFR activating mutations was observed in patients with progressive disease (PD, n=3).

Data presented here suggest that monitoring urine and plasma ctDNA for changes in levels of EGFR activating mutations during therapy may reflect tumor dynamics

Conclusions

 Mutation enrichment NGS testing by urine and plasma ctDNA correctly identified *EGFR* activating mutations in 85% of patients.

 Monitoring EGFR levels in urine/plasma enabled accurate assessment of response in advance of radiographic evaluation and regardless of therapy type in 100% of patients, with a cut-off of 90% decrease in EGFR load, discriminating between patients with disease control (PR or SD) and patients with disease progression (PD).

• Analysis of urine and plasma for *EGFR* mutations may be a viable approach for therapeutic monitoring of patients with advanced cancers.